

Ovation® WGA System (6100-12)

Enter the following information to automatically calculate the volumes needed to prepare each reaction. Simply print this page and record the essential information for this experiment.

Operator's Name: _____ Date: _____

Ovation WGA System Lot No. _____ RNAClean® Bead Lot Number _____

SPIA Product Purification Kit Name and Lot No. _____ Number of Samples _____

Thermal Cycler Programming	
PROGRAMMING DETAILS	
Program 1: Denaturation	95°C for 3 min 4°C forever
Program 2: Template Synthesis	4°C for 1 min, [23°C for 10 min, 35°C for 10 min, 57°C for 10 min] X 2, 70°C for 5 min 4°C forever
Program 3: SPIA Amplification	4°C for 1 min 47°C for 90 min 95°C for 5 min 4°C forever

Denaturation and Template Synthesis				
Thaw the Denaturation and Template Synthesis Reagents (blue) and Nuclease Free Water (green) . Mix each reagent, spin and place on ice				
Add 2 µL of DA1 into 0.2 mL PCR tube and place on ice				
For each sample add 8 µL of gDNA (10-50 ng), flick tubes to mix and spin				
Place the tubes in a thermal cycler running Program 1				
Immediately remove tubes from the thermal cycler after 95°C step and place tubes on ice for 5 min to snap cool				
Prepare Template Synthesis Master Mix (10% overfill is included in the calculation)	No. of Samples	DA2	DA3	DA4
	1	2 µL	7 µL	1 µL
Spin Template Synthesis Master Mix and place on ice				
Add 10 µL of the Template Synthesis Master Mix to each tube, mix and spin				
Place the tubes in a thermal cycler running Program 2				

Purification of Template cDNA

Bring the RNAClean beads to room temperature and resuspend

Add 32 μ L of beads to each reaction and mix

Incubate at room temperature for 10 min

Place tubes on magnet for 5 min to completely clear the beads

To minimize bead loss, remove only 40 μ L of Binding Buffer before the first wash step

Wash the beads while still on magnet for 30 sec with 200 μ L of freshly prepared 70% ethanol

Repeat wash 2 more times

Dry beads completely, at least for 15-20 min

Proceed immediately to SPIA® amplification.

SPIA Amplification

Thaw the **SPIA Reagents (red)**. Vortex **DB1** and **DB2**. Invert **DB3** 5 times, spin and place on ice.

Prepare **SPIA Master Mix** (10% overfill is included in the calculation)

No. of Samples	DB2	DB1	DB3
1	20 μ L	10 μ L	10 μ L

Mix **SPIA Master Mix**, spin and place on ice

On ice, add 40 μ L of **SPIA Master Mix** to each template synthesis reaction tube containing the dried beads, mix

Place tubes in a thermal cycler running **Program 3**

Once the thermal cycler reaches 4°C, spin and place tubes on ice

Note: Complete the remainder of the protocol in the Post-Amplification area

Purification of Amplified SPIA Product

Refer to the user guide and follow your method of choice for purification. Document the protocol below for future reference.

If not using RNAClean Bead method for the final SPIA product cleanup, place tubes on magnet plate and let stand for five min before transferring the cleared supernatant containing SPIA product to a fresh tube. Discard the beads

Proceed immediately to purification step or store SPIA product at -20°C

Add Binding Buffer in volume of:	Spin at speed:	For a duration of:
----------------------------------	----------------	--------------------

Add Wash Buffer in volume of:	Spin at speed:	For a duration of:
-------------------------------	----------------	--------------------

Repeat for second Wash.

To elute sample use **Nuclease-Free Water D1** provided with the Ovation WGA System

Add Nuclease-Free Water D1 in volume of:	Spin at speed:	For a duration of:
---	----------------	--------------------