

Protocol

Immunofluorescence Staining of Cells for Flow Cytometry

OVERVIEW

The following is a general cell staining protocol for analysis by flow cytometry. It is recommended that a separate isotype control be tested in parallel as a negative control.

MATERIALS REQUIRED

- Flow Buffer
- Trypan blue
- 10% formalin solution
- PBS
- 15 ml conical tubes
- 12 x 75 mm round-bottom tubes

MATERIAL PREPARATION

- **Flow Buffer**
 - 99 ml PBS (1X)
 - 1 ml FBS
 - 0.1 ml Sodium Azide (100%)
 - If detecting **extracellular antigens**, the Flow Buffer is ready for use.
 - If detecting **intracellular antigens**, add Saponin to a final concentration of 0.1%.
- Store Flow Buffer at 4°C.

STAIN CELLS FOR FLOW CYTOMETRY

1. Isolate and dissociate the cells to a single cell suspension. Collect the cells in a 15 ml conical tube.
2. Remove a sample of the cell suspension and count live cells using trypan blue and a hemacytometer.
3. Centrifuge the cell suspension at 300 x g for 5 minutes at 4°C.
4. Aspirate the medium and flick the tube to disrupt the cell pellet.
5. If detecting **extracellular antigens**, resuspend the cell pellet in 3 ml of PBS.
6. If detecting **intracellular antigens**, resuspend the cell pellet in 2 ml 10% formalin solution and incubate for 15 min at room temperature.

Note: Formalin should be added slowly while agitating the tube to avoid cell clumping.

7. Centrifuge the cell suspension at 300 x g for 5 minutes at 4°C.
8. Aspirate the medium and flick the tube to disrupt the cell pellet.
9. Add enough Flow Buffer to bring the cell suspension to a concentration of 2×10^6 to 1×10^7 cells per ml. Keep the cells on ice.

Note: The Flow Buffer formulation depends on if the antigen is extra or intracellular. Refer to the preparation section to make sure the appropriate Flow Buffer is used.



Stemgent
 One Kendall Square, B6201
 Cambridge, MA 02139
 Technical Support: 1.877.783.6911
 www.stemgent.com

Protocol

Immunofluorescence Staining of Cells for Flow Cytometry

10. For each sample, add 100 μ l of the cell suspension to a 12 x 75 mm round-bottom tube.
11. Add the appropriate amount of primary antibody and/or isotype control to each sample.
12. Incubate on ice for 30 minutes to 1 hour, protecting the samples from light if using a conjugated antibody.
13. Add 4 ml of Flow Buffer.
14. Centrifuge at 300 x g for 5 minutes at 4°C.
15. Aspirate the supernatant and flick the tube to disrupt the cell pellet.
Note: If using a conjugated antibody, skip steps 16 through 20 and go directly to step 21. If using a purified antibody, continue with step 16.
16. Add the appropriate amount of secondary conjugated antibody to each sample.
17. Incubate on ice for 30 minutes to 1 hour, protecting the samples from light.
18. Add 4 ml of Flow Buffer.
19. Centrifuge at 300 x g for 5 minutes at 4°C.
20. Aspirate the supernatant and flick the tube to disrupt the cell pellet.
21. Add appropriate volume of Flow Buffer to each tube.
*Note: For this step, use Flow Buffer **without** 0.1% Saponin for both extra and intracellular antigen detection.*
22. Analyze the cells by flow cytometry within 4 hours.



Stemgent

One Kendall Square, B6201

Cambridge, MA 02139

Technical Support: 1.877.783.6911

www.stemgent.com